
CAT

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Homogenising System X1740

Content

1	General Information	5
2	Unpacking the Instrument	5
2.1	Scope of delivery	5
2.2	Homogenising tools and accessories (to be ordered separately):	6
3	Safety Information	8
3.1	Explanation of the safety information in the instruction manual	8
3.2	Explanation of the safety information on the device	8
3.3	General Safety Information	9
4	Set-up and Operation	10
4.1	Intended Use.....	10
4.2	Description Drive Unit X1740	10
4.3	Set-up	11
4.3.1	Assembly of support rods	11
4.3.2	Attaching the drive unit to the U-stand	11
4.4	Description of a homogenising tool	12
4.4.1	Inserting the homogenizing tool	12
4.5	Operation X1740.....	13
4.5.1	Connect the unit to the mains.	13
4.5.2	Description of the speed scale:	13
5	Exchanging the homogenizing shafts.....	14
5.1	Exchanging the generator	14
6	Assembling the DK 40 Flow-Through-Chamber.....	15
6.1	Description of DK 40 Flow-Through-Chamber	15
6.2	Mounting the Flow-Through-Chamber to the Drive Unit X1740	16
7	Operation of Flow Through-Chamber.....	19
7.1	Flow-Through Mode	19
7.2	Circulation Mode	19
8	Maintenance and Cleaning	19
8.1	Cleaning the drive unit X1740	19
8.2	Cleaning the homogenising tool.....	19
8.3	Rough Cleaning of Flow-Through-Chamber and Tubing	19
8.4	Cleaning the Inside of the Flow-Through-Chamber and Generator (= Rotor and Stator)..	20
8.5	Maintenance of drive unit X 1740	21
8.6	Maintenance homogenizing tools.....	21

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1 General Information

The CAT homogenising drive unit X1740 is designed in accordance with Safety Class 1 and built and tested in accordance with DIN EN 61010.

According to these regulations the unit is designed to meet the requirements for safe and correct operations. To maintain the proper safety and operational functions of the instrument the user should follow the instructions and safety guidelines in this manual.

2 Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport to be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

2.1 Scope of delivery

Please check that the package contains the following:

Type	Description	Part No.
X 1740	Homogenizing Drive Unit, 230V	60408-00 or
X 1740	Homogenizing Drive Unit, 115V 2 support rods Instruction Manual	60408-01

2.2 Homogenising tools and accessories (to be ordered separately):

Tools	Description	Part No.
T 6	Shaft, 6 mm diameter, 55 mm long, complete with generator	60420*
T 6	Shaft, 6 mm diameter, 120 mm long, complete with generator	60410*
T 10	Shaft, 10 mm diameter complete with generator	60421*
T 17	Shaft, 10 mm diameter, complete with 17 mm generator	60426*
T 20	Shaft, 20 mm diameter complete with generator	60422*
G 20	Shaft, 20 mm diameter complete with generator	60423*
T 30	Shaft, 20 mm diameter complete with 30 mm generator	60424*
G 30	Shaft, 20 mm diameter complete with 30 mm generator	60425*
T 40	Shaft, 20 mm diameter complete with 40 mm generator	60414*
G 40	Shaft, 20 mm diameter complete with 40 mm generator	60415*
DK 40	Flow-through chamber with external cooling for G40 generator	
U-stand incl. rod	platform overall dimensions: 692x371 mm, rod diameter: 27mm, length:1010mm	60494-00
Special clamp	clamp to attach the unit to a stand	60495-00
Universal key	Universal rotor wrench for shafts 6 - 30 mm	60470-00
Socket wrench	for all shafts	60471-00
Rotor wrench	for shaft T 10	60472-00

**Important Note:**

Read this instruction manual carefully before operating the instrument. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer at the following address

CAT

Ingenieurbüro CAT
M. Zipperer GmbH
Etzenbach 16
D-79219 Staufen
+49-7636-7803-0

**Important Note:**

Please store the instruction manual in a place easily accessible to every user.

**Important Note:**

Make sure that the unit is to be operated under the conditions corresponding to that as described under the section titled "Technical Data"

A EURO AC plug (DIN 49441 CEE 7/VII10/ 16 A 250 V) is standard on all instruments. For North America instruments feature standard US plugs (NEMA Pub.No.WDI1961 ASA C 73.1. 1961 page 8 15A 125V).

When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.

The instrument is earthed as supplied. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

3 Safety Information

3.1 Explanation of the safety information in the instruction manual

This symbol emphasizes the following notices, which have to be strictly observed by the user! Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

This symbol emphasizes prohibitions, which have to be strictly observed. Any disregard of the prohibitions may result in severe impairment of serviceability as well as severe impairment of the user.

This symbol emphasizes notices, which have to be observed by the user to assure a safe operation of the unit.

CAUTION:
Hazard of fire or explosion!

3.2 Explanation of the safety information on the device

CAUTION:
This symbol indicates to read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked. Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

3.3 General Safety Information

Attention:

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

Attention:

Do not use the unit with flammable or explosive substances. The instrument is not explosion-proof.

Attention:

Only instructed users may operate the instrument.

Important Note:

Make sure that all safety information of the instrument respectively on the device itself are clearly visible during operation

Attention:

If noise level exceeds 85 dB (A) at the work station operator must use ear plug and/or ear mufflers.

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

Ensure that the unit is firmly attached to a stand.

Wear ear protection should the noise level exceed 85dB(A).


Wear eye protection.

4 Set-up and Operation

4.1 Intended Use


The X1740 is a homogenizing system which in connection with a homogenizing tool produces emulsions and dispersions and is used in laboratories, universities and industries. For safe operation the unit has to be fixed to a stand.

4.2 Description Drive Unit X1740

4.3 Set-up

4.3.1 Assembly of support rods


Screw the support rods (4) into the threaded bores (4a) at the back of the unit.

Insert a screw driver into the bores at the end of the support rods and firmly tighten the support rods by turning the screw driver to the right.

4.3.2 Attaching the drive unit to the U-stand

To ensure safe operation the drive units are mounted to the U-stand using the special clamp (8). Vibrations may loosen the fastening screws. Please check periodically whether all fastening screws are securely tight.

- Attach special clamp (8) to the stand rod.
- Now attach the drive unit to the stand by inserting the support rods (6) into the free openings of the special clamp (8).
- Tighten clamp screws (9) of the special clamp (8).


4.4 Description of a homogenising tool

Attention


The generators have sharp edges. Handle with care. Risk of injuries!

The homogenizing tool comes always assembled. The connecting adapter (drive pin (10)) is located on the top of the shaft and is used to attach the shaft to the drive motor.

4.4.1 Inserting the homogenizing tool

When inserting the homogenising tool, make sure that the coupling fits well. This can be easily detected by the "click" that is heard when the coupling shaft has reached the end of the coupling slot of the motor. Securely fasten the fastening screw (7).

Never run homogenising tools

dry as the bearings and gaskets will be damaged if the generator is not cooled by the medium. The distance between the homogenising tool and the vessel bottom should not be less than 10mm. The filling of the medium may not be less than about 55 mm. To avoid a vortex it is necessary to insert the dispersing shaft out of the center of the vessel. If these conditions are met the unit is ready for operation.

4.5 Operation X1740

Attention:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

4.5.1 Connect the unit to the mains.

The drive unit is turned on and off by means of the dip switch (1) at the motor housing.

- Dip switch position I: motor runs
- Dip switch position 0: motor off

With the speed adjustment wheel (3) the speed can be steplessly adjusted.

Caution:

The dip switch (1) is only single-pole and is thus not suitable for safely isolating the unit from the mains supply!

Before starting the unit make sure that

- the speed is set to the lowest speed
- the generator of the homogenising tool is inserted in the fluid. To avoid a vortex, it is necessary to insert the dispersing shaft out of the centre of the vessel.
- the vessel is protected against rotating, shaking or moving
- make sure that the air ventilation slots are always open, as this is necessary for cooling purposes

Start the unit always at the lowest speed and then gradually increase the speed to the desired level.


4.5.2 Description of the speed scale:

1	=	2500 rpm
2	=	7000 rpm
3	=	11000 rpm
4	=	15000 rpm
5	=	19000 rpm
6	=	23000 rpm

5 Exchanging the homogenizing shafts

Prior to exchange the homogenizing shaft disconnect the unit from the mains supply. Loosen the fastening screw (7) and firmly pull out the shaft by hand. It will disconnect easily. When inserting another shaft, make sure that the coupling fits well. This can be easily detected by the "click" that is heard when the coupling shaft has reached the end of the coupling slot of the motor. Do not forget to securely fasten the fastening screw (7).

5.1 Exchanging the generator

The generators have sharp edges. Handle with care. Risk of injuries!

To exchange the generator (rotor and stator = generator)(11) always take the shaft out of the motor drive unit as described above. For undoing the rotor, use the two special keys (to be ordered separately). The socket wrench is used to counter-hold the shaft, and the rotor wrench to undo the rotor. To exchange the stator, remember that the fastening thread is counter-clockwise. You will find further detailed descriptions in the instruction manual which comes with each homogenising tool.



Important Note

Always unscrew the rotor first before unscrewing the stator.


6 Assembling the DK 40 Flow-Through-Chamber

6.1 Description of DK 40 Flow-Through-Chamber


6.2 Mounting the Flow-Through-Chamber to the Drive Unit X1740


Take care to avoid damage of parts during assembling respectively disassembling.


Screw in the 3 thread bolts (19).


Attach the O-ring (20). Take care to uniformly press the O-ring (20) to flange DK40 (13).


Now insert shaft G20 (21)


Secure the shaft with threaded bolt (22)


Attach the chamber lid (23)


Turn the cooling jacket (18) and attach it to the flange of DK 40 (13)


Insert O-ring (24) and then screw the stator (25) to the shaft (21).
Attention: **left-hand thread!**

Now attach rotor (26).

Attach inlet chamber lid (17)

Secure inlet chamber lid (17)
with the 3 wing nuts (16).

Push motor flange (X1740) into the flange of DK40 (13). Turn the motor slightly back and forth until you feel that the motor is properly attached to the flow through chamber.

Note: The drive unit is firmly attached as soon as the wider part of the motor flange fits tightly to the flange of the DK40 (13). Turn the drive unit that way that the cooling slots at the motor flange show downwards. Now you can fix drive unit and flow-through-chamber with the help of the fastening screw (7).

7 Operation of Flow Through-Chamber

7.1 Flow-Through Mode

Connect the vessel containing the media to the inlet by a tube (inner diameter = 20 mm). Connect another tube (inner diameter = 20 mm) to the outlet and lead it to a collecting basin.

Note: Secure all tube connections with suitable tube clips.

7.2 Circulation Mode

Connect in and outlet tube as described in 6.1 and lead the outlet tube back into the original vessel. This way the medium might flow several times through the in-line unit.

Note: The flow-through-chamber is not self-priming. Inlet tube and flow through chamber have always to be filled with a medium. We therefore recommend to set up the vessel containing the sample at the same height as the flow through chamber and to additionally fit a blockage (for example a ball valve) to the connecting tube.

8 Maintenance and Cleaning

8.1 Cleaning the drive unit X1740

To clean the drive unit use only water with a detergent that contains tensides or use isopropyl alcohol for stubborn soiling.

8.2 Cleaning the homogenising tool

The generators have sharp edges. Handle with care. Risk of injuries!

To avoid clogging clean shaft, generator and in case of G-shaft the sealing parts after each use. This is done by operating it in a solvent which dissolves substance residues and is not harmful to the gasket. This is usually sufficient to clean the generator.

Chemical sterilization may be also a method. General-purpose disinfectants such as formalin, alcohol, etc. may be used. It is important to remove disinfectant residues with sterilized water.

Important Note:

Make sure that the bearings, O-rings and gaskets are resistant to solvents.

Sterilization with moist heat:

This denotes the use of a steam jet pressurized to 2 bars at 120°C.

8.3 Rough Cleaning of Flow-Through-Chamber and Tubing

Flow-through-chamber and tubing might be roughly cleaned by pumping a suitable liquid through the chamber instead of the medium.

8.4 Cleaning the Inside of the Flow-Through-Chamber and Generator (= Rotor and Stator)

- Disconnect first drive unit X1740 from the mains
- Detach the flow through chamber from the drive unit
- The flow-through chamber can be disassembled without disconnecting the tubes.
- Remove wing nuts (16) together with the washers.
- Now remove inlet lid (17). **Attention:** Rest liquid escapes!!!
- To loosen the rotor (26) you need the rotor wrench for shaft G40 plus the socket wrench for shaft G40.
Note: The wrenches have to be ordered as accessories.
- Push the socket wrench over the drive pin (10) and hold the socket wrench in this position.
- Insert the rotor wrench in the rotor and unscrew the rotor by turning to the left.
- Loosen the stator (25) by turning to the right and take it out. **Attention: Left hand thread.**
- The outer part of the slip ring sealing can now be taken out and cleaned.
- Please note the right sequence: Rotating slip ring, O-Ring, pressing disc and pressing spring.
- Now pull off the pipe of the cooling chamber with the chamber tube outlet (14) from chamber lid (23).

To re-assemble the DK40 please follow the steps described under Point 6

Please take care that neither O-rings nor any other parts experience any damage during assembling respectively disassembling the flow-through chamber.

Attention: Before operation check the slip ring sealing. The slip ring should be easily turned by hand. If it is stuck after a long period without use just loosen slip ring by wetting it with water. Never operate the flow-through chamber if the slip ring is stuck!!.

Note: Always unscrew the rotor before the stator.

Note: Rotor and stator wrenches are listed in the price list under accessories.

8.5 Maintenance of drive unit X 1740

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

8.6 Maintenance homogenizing tools

The generators have sharp edges. Handle with care. Risk of injuries!

The gaskets in the homogenizing tools must be constantly monitored. In the event of leakage the suction effect of the rotating shaft can cause the medium to penetrate as far as the drive unit. If liquid emerges from the side hole at the top of the shaft tube stop work immediately and check the gaskets. The function of the generators depends on the condition of the sharp edges on the rotor and stator. These edges may be blunted very quickly in abrasive media, reducing the effectiveness of homogenizing.

9 Dismantling, Transport and Storage

9.1 Dismantling

1. Switch the unit off.
2. Disconnect it from the mains supply.
3. Remove any glass beakers and any other equipment around the drive unit.
4. Remove the homogenizing tool by loosening the fastening screw (7).
5. Loosen the special clamp and remove the clamp from the support rod.
6. Unscrew the support rods.
7. Now you may remove the instrument from the working area.

9.2 Transport and Storage

Prior to transport:

Switch the instrument off and proceed with dismantling as described under "Dismantling".

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)
- Storage ambient:
 - Max. ambient temperature : RT to +40°C
 - Max. humidity: 80%

10 Disposal

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

11 Warranty and Liability

The manufacturer agrees to either repair, or replace, at the manufacturer's discretion, any defects in materials or workmanship which develop within 24 months of the delivery of this product to the original user. In the event of replacement, the replacement unit will be guaranteed for the remainder of the original twenty-four (24) months period or ninety (90) days, whichever is longer.

If this product should require service, contact your local distributor or manufacturer for necessary instructions.

This guarantee will not apply if the defect or malfunction was caused by accident, neglect, unreasonable use or fitness for a particular purpose, which extend beyond the description and period set forth herein.

The manufacturer's sole obligation under this guarantee is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind, resulting from use or possession of the product.

Attention: The user has to determine, if the instrument is suitable for his specific application. If there are any further queries, contact your local dealer or the manufacturer direct.

12 Technical Data

12.1 Technical Data X1740

Type: X 1740	Specifications
Motor Power input	1800 Watt
Motor Power output	1200
Idle speed	2500 - 23000 rpm
Voltage	230V/50 Hz or 110V/60 Hz
Case dimensions	100x100 x 360 mm, WxDxH
Weight	4.2 kg
Permissible ambient temperature	5 - 40 °C
Permissible humidity	80 % RH

13 Repairs

When returning instruments for repair that have come into contact with hazardous substances, please:

- Fill in attached "Repair Return Form"
- Provide precise information on the relevant medium
- Take protective measures to ensure the safety of our receiving and maintenance personnel
- Mark the package as appropriate for hazardous materials.

13.1 Repair Return Form

CONTACT/USER INFORMATION

Contact:	Phone No.
Fax No.	Email:
Billing:	Shipping:
Company	Company
Address	Address

INSTRUMENT INFORMATION

Model	Serial No.	
Please describe all problems/malfunctions		
Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to CAT and an MSDS for hazardous material must be included with the instrument.

RETURN SHIPPING

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	--	------------------------------------	----------------------------------

**Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

PACKAGING INSTRUCTIONS TO RETURN A CAT INSTRUMENT FOR REPAIR

- ✓ Remove all accessories (e.g. homogniser tools, stirring paddles) from the instrument
- ✓ Clean excess testing material off the instrument/accessory
- ✓ Include MSDS sheets for all hazardous materials used with this instrument
- ✓ Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- ✓ DO NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the instrument. (Applicable for overhead stirrers and homogenisers.)
- ✓ Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.

Example: "Handle with care" or "Fragile- Delicate Instrument"

14 Declaration of conformity

We declare under our sole responsibility that this product (see label) corresponds to the regulations and conforms to the standards or standardized documents

2006/95/EC

DIN EN 61010-1; VDE 0411-1:2002-08; Amendment 1; VDE 0411-1 Amendment 1:2002-11; Amendment 2; VDE 0411-1 Amendment 2:2004-01
DIN EN 61010-2-051; VDE 0411-2-051:2004-07

2004/108/EC

DIN EN 61326-1; VDE 0843-20-1:2006-10
DIN EN 61000-6-2; VDE 0839-6-2:2006-03
DIN EN 61000-6-4; VDE 0839-6-4:2007-09
DIN EN 61000-4-2; VDE 0847-4-2:2009-12
DIN EN 61000-4-3; VDE 0847-4-3:2008-06
DIN EN 61000-4-4; VDE 0847-4-4:2005-07
DIN EN 61000-4-5; VDE 0847-4-5:2007-06
DIN EN 61000-4-6; VDE 0847-4-6:2009-12
DIN EN 55011; VDE 0875-11:2010-05
DIN EN 55014-1; VDE 0875-14-1:2010-02
DIN EN 55014-2; VDE 0875-14-2:2009-06

2006/42/EC, 93/68/EEC

DIN EN ISO 12100-1:2004-04
DIN EN ISO 12100-2:2004-04

In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

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