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Table of Contents: Page 2

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Table of Contents

1. Scope of Delivery	5
2. General Information	5
2.1 Unpacking the Instrument	5
3. Safety Instructions	7
3.1 Description of Safety Symbols on the Device	7
3.2 Description of Safety Symbols in the Instruction Manual	7
3.3 Safety Instructions	8
4. Intended Use	9
5. Setting up the Instrument	9
5.1 Setup	9
5.2 Rear Connections	10
5.3 Mains Connection	10
5.4 Connection of the Pt100 Probe	10
5.5 Probe Stand	11
5.6 Extension Connectors	11
6. The Control Panel	12
6.1 Description of the Control Panel Functions	12
7. Heating	14
7.1 Temperature Control with external Pt100 Probe	14
7.2 Temperature Control without Pt100 Probe (Hot Plate Control)	15
7.3 Plate Limit (Heating with PT100-Probe)	15
7.4 Safety Temperature	16
8. Stirring	17
9. Fahrenheit-Celsius Setting	17
10. Timer Controlled Operation	18
11. Additional Safety Functions	18
11.1 Differential Alarm	18
11.2 Out of Liquid Check	18
11.3 Failure of External Pt100 Probe	18
11.4 Internal Over Heating Monitoring	19
12. Error Messages	19
12.1 Other Messages	20
13. Maintenance and Cleaning	20
14. Warranty and Liability	20
15. Technical Data	21

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Instruction Manual

Heatable Flask Stirrer KM16.4D / KM16.7D

- ✓ 4 Digit LED Display
- ✓ Possible Pt100 Temperature Control
 - ✓ Display of Set and Actual Values
 - ✓ User Friendly Control Panel
- ✓ PID Algorithms for Accurate Temperature and Motor Control
- ✓ Minimal Overshoot and High Temperature Stability
 - ✓ Motor Soft-Start
 - ✓ Programmable Heating Plate Limit
 - ✓ Programmable Safety Temperature
- ✓ Programmable Liquid Temperature in Pt100 Mode
 - ✓ Timer Shut-off Function
 - ✓ Over-Temperature Protection
- ✓ Shut Off Through Several Error Detection Mechanisms

1. Scope of Delivery

Unpack the delivered contents carefully and check to see that nothing appears physically damaged or is missing. Your delivery should contain the following:

- 1 Heatable flask stirrer KM16.4D (230V CAT Order No.: 60357-00, 110V CAT Order No.: 60357-01) or
- 1 Heatable flask stirrer KM16.7D (230V CAT Order No.: 60358-00, 110V CAT Order No.: 60358-01)
- 1 Instruction Manual
- 1 Set of centering rings for KM16.4D consisting of 3 parts or
- 1 Set of centering rings for KM16.7D consisting of 2 parts

Optional:

- 1 Pt100 temperature probe (CAT Order No.: 60278-00)
- Stand for KM16.4D (CAT Order No.: 20669-01)
- Stand for KM16.7D (only on demand)
- Probe Clamp (CAT Order No.: 8B005624)
- Clamp (CAT Order No.: 60668-00)
- RS485 cable for connection to PC (only on demand)

2. General Information

Our heatable flask stirrers are designed according to protection class 1. They were manufactured and tested according to DIN EN 61010.

According to these regulations, the units are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

2.1 Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

If the instrument is not damaged and all parts are complete you may start to operate the device after reading the instruction manual.

**Important Note:**

Read this instruction manual carefully before operating the instrument. Should there be any further questions after reading these instructions concerning the setup, operation or warranty please contact either your distributor or the manufacturer at the following address:

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**Important Note:**

Please store the instruction manual in a place easily accessible to every user.

A EURO AC plug (DIN 49441 CEE 7/VII10/ 16 A 250 V) is normally shipped with our instruments. North America instruments are shipped with standard US plugs (NEMA Pub.No.WDI1961 ASA C 73.1. 1961 page 8 15A 125V).

When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation. Ensure that your local supply voltage matches the indication on the instrument's specification plate!

The instrument is earthed. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

3. Safety Instructions

3.1 Description of Safety Symbols on the Device

Warning against general danger:

This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

Warning against hot surfaces: This symbol indicates hot surfaces on the device. It refers to the hotplate of the heatable flask stirrer which can reach up to 500°C during operation causing a hazard resulting in serious burns. Avoid contact with the hotplate, vessel and heating medium. Also avoid storage of flammable substances near the instrument.

Caution!

Beware of the effect of the magnetic field on cardiac pacemakers and data storage media.

3.2 Description of Safety Symbols in the Instruction Manual

This symbol points to the following notices, which have to be strictly observed by the user! Disregard of warnings may result in the impairment of serviceability as well as in physical harm to the user.

Caution

Danger for fire or explosion!

Warning against hot surfaces:

This symbol indicates that there are hot surfaces on the device. It refers to the hotplate of the magnetic stirrer which can reach up to 500°C during operation causing a hazard resulting in serious burns. Avoid contact with the hotplate, vessel and heating medium. Also avoid storage of flammable substances near the instrument.

Note regarding repair or maintenance

Note regarding mains cable connection

Caution

Note regarding mains voltage

This symbol emphasizes notices which have to be observed by the user to ensure safe operation of the unit.

3.3 Safety Instructions

Caution:

Please comply with all safety and accident-prevention regulations applicable to laboratory work!

Caution:

Use great caution when working in the vicinity of flammable or explosive substances. This instrument is not explosion-proof.

Caution

Instructed users only may operate the instrument.

Caution:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.

Note:

Make sure that all safety information on the instrument is clearly visible during operation.

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

Connect the instrument to an earthed AC power outlet only.

Attention:

When using flammable substances set the temperature to at least 25°C lower than the flashpoint of the sample.

Ensure that the unit is standing on a solid and horizontal surface.

4. Intended Use

These heatable flask stirrers are mainly to be used to heat and stir hydrous solutions and oil not containing inflammable parts in glass vessels. They are designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

5. Setting up the Instrument

5.1 Setup

Please place the unit on a fire-proof horizontal and even surface. A minimum distance of 50 cm to inflammable materials should always be observed.

**Caution**

The power cable may not come in contact with the hotplate.

**Caution**

The unit may not be operated in explosion-prone areas.

**Caution**

The unit is not to be used without supervision.


**Caution**

When connecting the instrument to an AC power outlet ensure that your local supply voltage matches the indication on the instrument.

**Caution**

Beware of the effect of the magnetic field on cardiac pacemakers and data storage media.

5.2 Rear Connections

(Rear View of Magnetic Stirrer)


5.3 Mains Connection

Connect the mains Power Cable to the mains power connector. Ensure that your local supply voltage matches the indication on the instrument

5.4 Connection of the Pt100 Probe

A DIN Connector on the rear of the instrument is intended for the connection of a Pt100 temperature Probe (CAT Order No. 60278-00). Pt1000 Probes are not suitable.

Verbindung eines PT100 Sensors:

Note:

Use Pt100 probes with electrically isolated tips only. Probes that allow for an electrical connection to the metal chassis of the unit may distort temperature measurement.

Ensure that the Pt100 probe tip is submersed to a depth of at least 50mm in the heating medium to allow for accurate measurement.


Caution

Ensure that the cable of the Pt100 probe does not come in contact with the heating plate.

5.5 Probe Stand

Use the M6 Thread on the unit for the optional probe stand (CAT Order Nr. 20669-01).

5.6 Extension Connectors


These are intended for the connection of external modules or RS485-RS232 adapters for RS485 communication to a Personal Computer. Please contact your supplier or manufacturer for more information and order numbers.

Please note that these connectors are intended for CAT accessories and cables only. Connection of other cables (e.g. Network cable) and devices will cause permanent damage to the device and to other connected electronics.

6. The Control Panel

The control wheel on the left is used for changing set values. Any values and messages to be displayed are shown on the 4 digit 7-segment display. The turning knob on the right is used for setting the speed of the stirring motor. When switching on the unit the set motor speed will be displayed on the LED.

(Control Panel of the Magnetic Stirrer)

6.1 Description of the Control Panel Functions

7Segment LED Display

This LED display is used for showing different values and messages (temperature, motor speed, timer and error messages)

Standby/On Button

Switch the device on or off (standby). Switching the unit off all set temperature values will be stored internally.

Button: „Plate On/Off“

To switch on the heating plate after entering the desired set-values this button needs to be pressed. Switching the plate off is also done by pressing this button.

Button: „View RPM“

Show the actual set motor speed

**Button: „Plate Temp“**

Display the actual plate temperature or press again for setting the set point temperature. When used with Pt100 probe: display the actual plate temperature or press again for setting the plate limit temperature.

**Button: „Safety Temp“**

Display or set the safety temperature

**Button: „Probe Temp“**


Display or set the probe temperature (Pt-100)

**Button: „Timer“**

Display or set the remaining timer minutes

**Motor Speed Wheel**

This wheel is used to set the speed of the stirring motor. The set speed is displayed on the LED display while this wheel is turned (opto-electronic speed control)

**Encoder Wheel with Push Button**

This encoder wheel is used to change different set values. Turning clockwise will increment values, while turning anti-clockwise will decrement. Pressing the encoder wheel push button will confirm a set value. Any error messages displayed at start can also be reset by pressing the encoder wheel push button.

Status LED „Plate On/Off“

This LED is lit when the hotplate is on. It indicates that the unit is trying to reach or keep the entered set temperature for hotplate or probe.

Status LED „Heater“

This LED indicates that the hotplate is currently being exposed to an energy pulse.

Status LED „Hot“

When the hotplate is hotter than 60°C this LED lights up.

Status LED „Timer“

When the timer is activated this LED lights up.

7. Heating

The heatable flask stirrer allows for 2 different modes of operation:

- **Heat control with external Pt100 probe temperature control (Recommended)**
- **Heat control with hotplate only (no Pt100 probe connected)**

For exact temperature control it is recommended to use an external Pt100 probe.

In order to avoid the round bottom flask falling over you should use the matching centering ring.

Caution: The device will remain hot after switching off the device!

If the power cable is connected and the unit is turned off (standby) the display will show the message „HOT“ as long as the hot plate temperature exceeds 60°C.

!!!Always switch off the device after use!!!

7.1 Temperature Control with external Pt100 Probe

We recommend the use of an external Pt100 probe (CAT Order Nr. 60278-00) to allow for accurate temperature control of the heating medium.

- **Connect** the Pt100 probe with the DIN connector on the rear of the heatable flask stirrer unit.
- **Place** the Pt100 probe in the vessel (min. 5cm submersion)

Press Power On/Stby to switch on the device. Any error messages displayed at start can be reset by pressing the encoder wheel button.

When „**A**“ (Amount) is shown as the first digit in the display **turn the encoder wheel** until the desired volume of the heating medium is shown. Press the encoder wheel button to confirm the volume in the vessel.

Press Probe Temp (press again if the **P** does not blink). Set the desired setpoint temperature value by turning the encoder wheel. Press the *Probe Temp* key again or press the encoder wheel button to confirm the set value or wait 3 seconds.

Press Plate On/Off to switch on the hotplate (Status LED above the button will now be lit).

When a Pt100 probe is connected, the value set with „Plate Temp“ does not refer to setpoint of the hotplate but to the plate limit temperature (Plate Limit – See section 7.3).

Normally when using a Pt100 probe the plate limit should be set to its maximum value (KM16.4D and KM16.7D: 450°C) to ensure optimal temperature control.

When the Pt100 probe is connected to the device, the setpoint temperature of the probe will be set automatically to 0°C, the hotplate will be switched off and the plate limit (see section 7.3) temperature will be set to its maximum value.

Caution

Please ensure that the cable of the Pt100 probe does not come in contact with the hotplate.

7.2 Temperature Control without Pt100 Probe (Hot Plate Control)

For certain applications it could be necessary to control the hotplate temperature instead the probe temperature. This is only possible if no Pt100 probe is connected to the device.

Press **Power On/Standby** to switch on the device. Any error messages displayed at start can be reset by pressing the encoder wheel button.

Press **Plate Temp** (press again if the **H** does not blink). Set the desired setpoint temperature value of the hotplate by turning the encoder wheel. Press the **Plate Temp** key again or press the encoder wheel button or wait 3 seconds to confirm the set value.

Press **Plate On/Off** to switch on the hotplate (Status LED above the button will now be lit).

If a Pt100 probe is connected to the device the set value of the “plate temp” is not the setpoint of the hotplate temperature but the maximum plate temperature limitation.

After disconnecting the Pt100 probe the setpoint of the hotplate temperature will be set automatically to 0°C.

7.3 Plate Limit (Heating with Pt100 probe)

For normal operation, the plate limit temperature should be set to the maximum value to achieve optimal heat-up times. Setting the limit temperature to the maximum value allows the microprocessor to freely select the optimum hotplate temperature to reach the set probe-temperature. However, when operating with flammable liquids it may be required to limit the maximum allowed hotplate temperature to a safe value.

Note: If the hotplate temperature is limited to a too low value the set temperature for the fluid cannot be reached or an extremely long heat-up time will be the case.

If the setpoint of the external Pt100 probe is changed, the plate limit temperature will automatically be set to its maximum value. The new value appears in the display and could be changed to the desired value directly by using the encoder wheel.

To set the hotplate limit temperature when Pt100 probe is connected:

Press **Plate Temp** (press again if the **H** does not blink). Set the desired plate limit temperature by turning the encoder wheel. Press the *Plate Temp* key again or press the encoder wheel button or wait 3 seconds to confirm the set value.

7.4 Safety Temperature

A safety temperature can be set to shut down the device if there is an error or to prevent an unexpected heat up of the medium.

The safety temperature is set automatically by the unit when the Pt100 probe or hotplate setpoint temperature is changed. If the unit is operated with an external Pt100 probe and the setpoint is changed the device automatically set the safety temperature 15°C above the new setpoint. If the liquid temperature exceeds the safety temperature (temperature measured with the external probe) the unit will shut down.

Conversely, if the unit is operated without an external probe and the hotplate setpoint is changed the device automatically set the safety temperature 15°C above the new hotplate setpoint. If the temperature measured at the hotplate (internal sensor) exceeds the safety temperature, the unit will also shut down.

If the device automatically changes the safety temperature the new value is shown in the display for a short time. After the device automatically set the new value another safety temperature can be set manually by user.

To set the safety temperature:

Press **Safety Temp** (the **L** starts to blink). Set the desired Safety temperature by turning the encoder wheel. Press Safety Temp again to confirm the set value or press the encoder wheel button or wait 3 seconds. The set value will be effective until setpoint of *Probe Temp* or setpoint of *Plate Temp* are changed and the *Safety Temp* is automatically set.

When a Pt100 probe is used the safety temperature refers to the temperature measured by the Pt100 probe. Otherwise the safety temperature refers to the hotplate.


In safety critical applications the proper use of the safety temperature is advised. Whenever the setpoint temperature of the external probe (Pt100) or the hotplate is changed the current safety temperature is set to 15°C above the setpoint of the *Probe Temp* or the *Plate Temp* value. To set the safety temp as close to the actual setpoint value as possible, set the probe or plate temperature first, before changing the safety temperature value.

8. Stirring

Make sure that the vessel used is permeable to magnetic field lines (e.g. glass, ceramics or stainless steel)

To set the motor speed:

Set the motor speed by turning the motor speed wheel

The set motor speed is shown on the LED display while the motor speed wheel is being turned. Turn the wheel entirely to the left to switch the motor off.

Note: Optimal magnetic coupling is acquired with stirring bars with lengths: 30-60mm. Use stirring bars that are suitable for magnetic flask stirrers (egg shaped).

Caution

Use a suitable vessel (round flask) and beware of high stirring speeds and high liquid levels to avoid splashing. Make sure that suitable protective clothing and eyewear is used.

9. Fahrenheit-Celsius Setting

This heatable flask stirrer can show temperature values in degrees Celsius as well as Fahrenheit. When the device is switched on in Celsius mode the normal start sequence will be displayed. In Fahrenheit Mode the message **FAHr** will be displayed after switching on the device.


You can change between Celsius and Fahrenheit mode in the following way:

Press **Power On/Stby** to switch on the device.

Shortly **after** switching the device on **press the Probe Temp Key and hold**.

The display now shows a flashing **U** and you can change between the settings **CEL** (Celsius) or **FAH** (Fahrenheit) by turning the encoder wheel. Press the encoder wheel to confirm the selected temperature mode. The unit will remember the selected temperature mode until it is changed again by the user.

10. Timer Controlled Operation

This device has an internal timer function. The timer allows for automatic switching off. The switch off time can be selected in the following way.

E.g. a set time of 10 minutes will cause the shut-off of the unit after 10 minutes.

Press **Timer** key (press again if the  does not blink). Turn the **encoder wheel** to select the desired switch off time (remaining on-time in minutes). Press the Timer key again or press the encoder wheel button or wait 3 seconds to confirm.

Minimum time until Switch off:	1 minute
Maximum time until Switch off:	999 minutes = 16h 30min

11. Additional Safety Functions

The following safety functions provide additional protection from uncontrolled heating. If a safety problem is detected, the device turns off automatically. The reason for the shutdown is displayed on the LED-display. When switching on the device next time, you have to confirm the displayed safety problem by pressing the encoder wheel button.

11.1 Differential Alarm

Breakage of the vessel could create a hazardous situation. To monitor this situation problem this instrument has a '*differential alarm*' feature. The '*differential alarm*' switches the unit off if the probe has fallen out of the liquid (e.g. breakage of glass). The '*differential alarm*' operates if the temperature fall detected by the Pt100 probe exceeds the internal differential alarm sensitivity. The unit then switches off.

11.2 Out of Liquid Check

The unit monitors, whether the Pt100 probe tip is immersed in the liquid in the vessel on the hotplate. The microprocessor checks whether the probe temperature changes in relation to the temperature changes of the hotplate. If the microprocessor detects no increase in the liquid temperature over a certain time, although the hotplate temperature is rising, a warning message is shown on the display and the unit shuts down.

Note:

This safety feature is only active if the setpoint temperature of the probe is at least 10°C higher than the actual measured temperature.

11.3 Failure of External Pt100 Probe

Disconnection or failure of the external Pt100 probe after the device is switched on will be detected and cause the unit to switch off.

11.4 Internal Over Heating Monitoring

The temperature inside the device is constantly monitored. A safety mechanism will switch the device off in case the internal temperature exceeds a safe value.

12. Error Messages

This unit has an on-board self-diagnostic program that detects possible failures. In case of an error message the program switches off the unit. To enable the user to prevent the problem from re-occurring, the unit will display any of the following error messages as it is switched on again. To reset any error message press the encoder wheel button:

LED Display:	Error name:	Description/Reasons :
HErr	hotplate failure	The hotplate temperature exceeded its maximum value The hotplate temperature sensor may be broken
OFEr	Off Error	The device was not switched off with the power button. Disconnecting the power cable or a mains power failure may be the source for error.
H4AF	Plate > Safety	The measured hot plate temperature was higher than the set safety temperature for the hotplate (Safety Temp.)
P4AF	Probe > Safety	The measured probe plate temperature was higher than the set safety temperature for the Pt100 probe (Safety Temp.)
CLOC	TIMER expired	Unit switched off after timer expired (→ no error, user information only).
dI FF	Differential Alarm	The differential alarm was triggered The temperature drop in the probe was too quick.
PFAI	Probe failure	The Pt100 probe is broken or damaged The probe was disconnected during operation.
t-r-I-A	Triac failure	The heating plate output stage is damaged.
HFAI	Plate Sensor Fail	The hotplate sensor is broken
HOU4	Housing overtemp	The internal temperature exceeded its maximum safe value
nLI 9	Out Of Liquid	The external temperature probe is not inserted into the medium.
COnt	Contact Thermometer Error	The contact thermometer is damaged or was disconnected during operation.
Udt0	Watch Dog Timer	The Watch Dog Timer was triggered.

12.1 Other Messages

LED Display:	Error name:	Description/Reasons :
FAhr	Fahrenheit	°Fahrenheit Mode
OFF	OFF	Unit was switched off by ON/OFF key
4850	RS232 Off	Unit was switched off by RS232/RS485 PC Communication
HOT	Hotplate is hot	The device is switched off but the hotplate exceeds 60°C.

13. Maintenance and Cleaning

The outer casing is mostly 304 grade stainless steel, the hotplate consist of aluminum. Together with a chemically resistant splash-proof membrane key-pad and therefore easily cleaned with warm water and any proprietary liquid laboratory detergent. Do not use steel wool or any similar plastic wool sponge or any aggressive cleaning agents to clean the unit.

Crust, limestone and water stains may be cleaned with the usual household detergents. Please remove remaining detergent thoroughly with a damp cloth or sponge, since some of the detergents react caustic with higher temperatures. If soils are already encrusted use a razorblade scratcher. Plastics, kitchen foil as well as stuff containing sugar must be removed immediately.

Mechanical supports are: Razorblade scratchers, non-scratching sponges and non-scratching steel wool.

The KM16.4D and KM16.7D should only be opened and repaired by authorized service personnel. Any work on the electronics in the unit should only be carried out by knowledgeable, trained personnel. Any attempt by the user to repair the unit will immediately render the guarantee null and void. Please contact your local distributor in the event of a problem.

In case of malfunction do not attempt to repair the unit. There are no user-serviceable parts in this instrument.

14. Warranty and Liability

The manufacturer agrees to correct for the original user of this product, either by repair, or at the manufacturer's election, by replacement, any defects in material or workmanship which develop within 24 months after delivery of this product to the original user. In the event of replacement, the replacement unit will be warranted for the remainder of the original twenty-four (24) month period of ninety (90) days, whichever is longer.

If this product should require service, contact your distributor or manufacturer for necessary instructions.

This warranty shall not apply if the defect or malfunction was caused by accident, neglect, unreasonable use, improper service, or other causes not arising out of defects in material or workmanship.

There are no warranties, expressed or implied, including, but not limited to, those of merchantability or fitness for a particular purpose, which extended beyond the description and period set forth herein. The manufacturer's sole obligation under this warranty is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damage of any kind resulting from use or possession of the product.

Note

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

15. Technical Data

Type KM16.4D, KM16.7D	Specifications
Mains	110/230 Volt; 50/60 Hz (See rear panel)
Display	4 digit 7Segment-LED display
Motor	<ul style="list-style-type: none"> - split-pole, 25 Watts - 60 - 1100 min⁻¹ - Softstart, microprocessor controlled
Temperature control	<ul style="list-style-type: none"> - microprocessor controlled (PID Control) - Pt100 probe adjustable in 1°C Steps - Hotplate adjustable in 1°C steps up to 450°C - programmable probe temperature (Pt100) - programmable hot plate temperature - programmable safety temperature
Safety system	detects and protects against the following hazardous situations: -hotplate failure (hotplate exceeds max. allowed temperature) -probe failure (disconnection or break of Pt100 probe) -detection of out-of-liquid condition, e.g. stirring vessel breaks -differential safety cut-out (external probe senses an extreme drop in temperature of liquid being stirred/heated) -triac failure
Switch-off time delay (Timer)	programmable switch-off time delay min. delay: 1 minute max. delay : 999 minutes
Temperature sensors	Hotplate: Thermocouple External: Platinum Resistance Sensor (Pt100)
Setpoint resolution hot plate	1°C / 1K (1°F in Fahrenheit Mode)
Setpoint resolution Pt100 Probe	1°C / 1K (1°F in Fahrenheit Mode)
Motor speed resolution	10 min ⁻¹

Typical temperature stability with Pt100 use*	± 0.4 °C [†]
Temperature range heating medium	Up to 330°C
Recommended stirring bar	Egg shaped
Permissible Ambient Temperature	5 - 40 °C
Permissible Humidity	80 % RH
Safety class acc to DIN 40050	IP32
Case dimensions	KM16.4: 200 x 145 x 175 mm KM16.7: 200 x 190 x 195 mm
Weight	KM16.4: 3.2 kg KM16.7: 3.4 kg

Type KM16.4D	
Hot plate material	Aluminum
Hot plate dimensions	Ø 90 mm round
Heating power	500 Watts

Type KM16.7D	
Hot plate material	Aluminum
Hot plate dimensions	Ø 120 mm round
Heating power	500 Watts

* 2L H₂O at 60°C Setpoint, Ambient: 23°C. Maximum variance during 60 min test.

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