**CAT**

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11.1	CHEMISCHE BESTÄNDIGKEIT:.....	26
<b>12</b>	<b>GARANTIE UND HAFTUNGSAUSSCHLUSS .....</b>	<b>27</b>
<b>13</b>	<b>ERKLÄRUNG ZUR GESUNDHEITLICHEN UNBEDENKLICHKEIT .....</b>	<b>28</b>
<b>14</b>	<b>INSTRUCTION MANUAL.....</b>	<b>29</b>
<b>15</b>	<b>GENERAL INFORMATION .....</b>	<b>30</b>
15.1	UNPACKING THE INSTRUMENT .....	30
15.2	EXPLANATION OF THE SAFETY INSTRUCTIONS IN THE OPERATING MANUAL .....	31
15.3	SAFETY INFORMATION.....	31
<b>16</b>	<b>PURPOSE OF THE CONTIBURETTE.....</b>	<b>33</b>
16.1	OPERATING EXCLUSIONS .....	34
<b>17</b>	<b>DESCRIPTION OF THE CONTIBURETTE .....</b>	<b>35</b>
17.1	OPERATING ELEMENTS .....	35
17.2	EXPLANATION OF FRONT PANEL PUSH LEVERS .....	35
17.3	DESCRIPTION OF THE HANDWHEEL.....	35
17.4	DISPLAY .....	36
<b>18</b>	<b>SET-UP .....</b>	<b>36</b>
18.1	INSTALL SUCTION TUBE.....	36
18.2	FITTING BURETTE TO BOTTLE .....	37
18.3	DISCHARGE TUBE ASSEMBLY .....	37
18.4	BATTERY INSTALLATION .....	37
<b>19</b>	<b>OPERATING OF THE CONTIBURETTE.....</b>	<b>38</b>
19.1	SWITCH ON/OFF .....	39
19.2	FILLING THE INSTRUMENT AND REMOVING THE AIR .....	40
19.3	OPERATING MODES OF THE CONTIBURETTE .....	40
19.3.1	<i>Titrate-Mode</i> .....	40
19.3.2	<i>Dosing-Modus</i> .....	41
19.4	SETUP MENU .....	41
19.4.1	<i>Setting of the display backlight</i> .....	41
19.4.2	<i>Setting of the automatic switching off</i> .....	42
<b>20</b>	<b>CALIBRATION OF THE CONTIBURETTE .....</b>	<b>43</b>
20.1	PROCEDURE OF CALIBRATION:.....	43
<b>21</b>	<b>ADJUSTING OF THE CONTIBURETTE.....</b>	<b>44</b>
21.1	PROCEDURE TO ADJUST THE CONTIBURETTE (USERCAL.): .....	44
21.2	RETURN TO FACTORY SETTING .....	45
<b>22</b>	<b>CLEANING AND MAINTENANCE .....</b>	<b>46</b>
22.1	PREPARATION FOR CLEANING .....	46
22.2	CLEANING PROCEDURE .....	47
22.3	CLEANING THE CONTIBURETTE.....	47
22.4	STERILIZATION .....	47
<b>23</b>	<b>DISMANTLING AND DISPOSAL.....</b>	<b>48</b>
23.1	DISMANTLING .....	48
23.2	DISPOSAL .....	48

<b>24</b>	<b>TRANSPORT AND STORAGE</b> .....	<b>48</b>
24.1	TRANSPORT/STORAGE .....	48
24.2	RETURN FOR REPAIR OR CALIBRATION.....	49
<b>25</b>	<b>WARRANTY AND LIABILITY</b> .....	<b>50</b>
<b>26</b>	<b>TECHNICAL DATA</b> .....	<b>51</b>
26.1	CHEMICAL RESISTANCE: .....	51
<b>27</b>	<b>DECLARATION ON THE ABSENCE OF HEALTH HAZARDS</b> .....	<b>52</b>

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**$\mu$  20 H**

## 15 General Information

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport be reported at the time of unpacking. Notify your supplier and the carrier or forwarding agent immediately in case of such damage.

### 15.1 Unpacking the instrument

Unpack the outer carton containing your instrument. If there are any punctures, side walls crushed, extensive chemical stains, water marks or other physical evidence that the contents may have been damaged, notify the carrier of the potential problem and ask for specific instructions. If there is no visible damage to the carton's exterior, open it in the normal manner. When all parts are laid out, you should have:

- 1 Contiburette  $\mu$  xx H (xx depends on type of the Contiburette)
- 1 Standard discharge tube
- 2 Adapters (28mm, and 45mm)
- 1 Filling/suction tube (May be cut to desired length)
- 2 1.5 Volts AA alkaline batteries
- 1 Instruction Manual

A selection of additional adapters, discharge tubes and accessories may be found in the pricelist.

**Note:**

Read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked.

Should you have any additional questions, after reading these instructions, concerning the setting up, operation or warranty terms, please contact either your distributor or the manufacturer at the following address:

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## 15.2 Explanation of the safety instructions in the operating manual

	This symbol highlights the following instructions, which the user has to follow strictly! Failure to follow these instructions can impair the safe function and safety of the user.
	This symbol highlights bans, which the user has to follow strictly! Failure to follow these bans can seriously impair the functionality and safety of the user.
	This symbol highlights instructions that should be strictly followed by the user to ensure safe operation of the device.

## 15.3 Safety Information

	<b>Attention:</b> Extreme caution should be exercised when handling toxic, corrosive, fuming, volatile or any other potentially dangerous substances.
	<b>Danger:</b> Do not use the device in potentially explosive areas. Mortal danger!
	<b>Attention:</b> Use great caution when working in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof. Do not use flammable or explosive substances near the instrument.
	<b>Attention:</b> Do not open the instrument. Repairs are only to be carried out by trained service technicians.
	<b>Attention:</b> Observe all markings on the reagent bottles. Dangerous and fuming chemicals must be dispensed in a fume hood. Only employ the instrument for the purpose intended by the manufacturer, and particularly within the resistance limits of the instrument. If in doubt, contact your supplier, or the manufacturer's factory representative at the phone number shown at the front page of this operating instruction.
	<b>Attention:</b> Always use the instrument in such a manner that neither the operator, nor any other person is endangered. When dispensing, maintain a physical distance between the instrument and the body. Avoid splashes; use proper connecting vessels, protective clothing and gloves.

	<p><b>Attention:</b> Please comply with all safety and accident-prevention regulations applicable to laboratory work.</p>
	<p><b>Attention:</b> Only instructed users may operate the instrument.</p>
	<p><b>Attention:</b> Always check the instrument for leaks and air bubbles. Special attention should be directed to determine that all push-ons, threaded connections and suction tubes are firmly in place before beginning operation. Leaking solutions may endanger persons and materials</p>
	<p><b>Note:</b> This instruction sheet does not purport to address all of the safety problems which might result from the use of this instrument, chemicals, reagents, apparatus or equipment employed in any specific test or protocols. It is the responsibility of the user to consult and establish appropriate safety and health practices, and then determine the applicability of regulatory limitations prior to use.</p>
	<p><b>Note:</b> Read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked.</p>
	<p><b>Note:</b> Please store the instruction manual in a place easily accessible to every user.</p>
	<p><b>Note:</b> Only use original spares and accessories. These are in alignment with the instrument and ensure reliability.</p>
	<p><b>Note:</b> In case of trouble (e.g., piston difficult to move or leakage), immediately stop pumping. Clean the instrument according to the cleaning instructions before any further use of the instrument or contact the manufacturer.</p>
	<p><b>Note:</b> Please make sure that all safety signs, also on the instrument, are visible.</p>

## 16 Purpose of the Contiburette

	<b>Attention:</b> Please comply with all safety and accident-prevention regulations applicable to laboratory work.
	<b>Attention:</b> Do not use flammable or explosive substances near the instrument.
	<b>Attention:</b> It is the responsibility of the user to consult and establish appropriate safety and health practices, and then determine the applicability of regulatory limitations prior to use. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer.
	<b>Attention:</b> Use the instrument only in compliance of the intended purpose and in way that neither user nor any other persons are endangered.

This instrument is designed for dosing liquids up to a concentration of max. 2 Mol/l, observing the following physical limits:

- 15 to 40 °C of instrument and reagent
- When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials:
  - Al<sub>2</sub>O<sub>3</sub>,
  - ETFE

During operation, the Contiburette must be screwed on a bottle, standing on a plain solid basis. It might be necessary to lock the bottle against tilting.

## 16.1 Operating Exclusions

	<b>Attention:</b> Compatibility of the instrument for the application must be checked by the user or contact the manufacturer
	<b>Attention:</b> Do not use flammable or explosive substances near the instrument.
	<b>Attention:</b> It is the responsibility of the user to consult and establish appropriate safety and health practices, and then determine the applicability of regulatory limitations prior to use. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer.


### Never use this instrument for

- liquids attacking Al<sub>2</sub>O<sub>3</sub>, ETFE
- suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument
- strongly crystallizing solutions, concentrated acids and bases as well as non-polar solvents which effect swelling of ETFE
- CARBON DISULPHIDE, AS THIS MEDIA INFLAMES EASILY

The Contiburette must not be autoclaved!

## 17 Description of the Contiburette

### 17.1 Operating elements

### 17.2 Explanation of front panel push levers

- On:** Push button to switch on the Contiburette
- Off:** Push button to switch off the Contiburette
- Left / - und Right / +:** Push button to change the menu and to modify a parameter. Permanent pressing effects a quick change of the value.
- Set Value / OK:** Button to enter a menu and to confirm a modified parameter
- Back / Reset:** Button to leave a menu and to refuse a changed value

### 17.3 Description of the handwheel

Turning the handwheel effects dosing of the liquid out of the bottle. Redosing is not possible. The dosed amount is only calculated, when the device is switched on.

## 17.4 Display

**cal.!** **Volume: 5,83 ml**  
**Dose: 10,00 ml**      **Rest: 4,17 ml**  
**Titrate** **Dose** **Setup** 

<b>Calibration status:</b>	In the left upper corner of the display the actual calibration status is shown. No indication means, that the device is factory calibrated. When user calibrated <b>cal.!</b> is displayed.
<b>Dosed Volume:</b>	The dosed volume is displayed in the right upper corner after <b>Volume:</b>
<b>Target volume and rest volume:</b>	The set volume (in the dosing mode) is shown after <b>Dose:</b> (When no volume is set you see: <b>Cont.</b> ) The remaining quantity is shown after <b>Rest:</b>
<b>Menu bar:</b>	The selected menu in the menu bar in the lower part of the display is accentuated with a black background.
<b>Batterie status:</b>	In the right lower corner you see the actual charge state of the batteries. When the capacity is lower than 10% <b>low Bat.</b> is displayed. The display backlight is disabled to save energie. The Contibürette is switched off automatically when the batteries are empty.

## 18 Set-up

Ensure that the instrument is standing on a solid surface.  
Please observe the specified ambient conditions (temperature and humidity) listed under „Technical Data“ as well as the safety instructions.

### 18.1 Install suction tube

**Note:**

**Loose connections such as incomplete push-ons, loose threaded fittings or a poorly fitting suction tube lead to ventilation in the system. Inaccurate measurements will result !**

Push an end of the filling tube onto the lower end of the pump housed in the base of the burette. It should be pushed as far as it will go. Then cut off the tubing at the bottom end to the length desired appropriate for the bottle size. The cut should be at an angle to prevent suction blockage.

## 18.2 Fitting burette to bottle

When fitting the burette to the titration bottle do not hold the burette on the black hand wheel - hold it on the grey housing instead. Do not carry the mounted instrument by its housing. Carry the mounted instrument by placing one hand underneath the bottle and the other hand around the base of the instrument at the bottleneck.

Two adapters have been provided with your instrument. These adapters enable the burette to fit into bottles of varying neck sizes. Select the appropriate adapter and twist the male threaded end snugly into the base of the burette. Other adapters are available.

## 18.3 Discharge tube assembly

**Attention:**

Always check the Contiburette for leaks and air bubbles. Special attention should be directed to determine that all push-ons, threaded connections and suction tubes are firmly in place before beginning operation.

**Leaking solutions may endanger persons and materials.**

Install the tube provided, by pushing it gently onto the threaded post at the left of the burette as far as it will go. Screw the rigid collar firmly into place and orient the tip of the discharge tube to the beaker, flask or other container receiving the dispensed liquid.

## 18.4 Battery installation

Two AA batteries provided are to be installed at the back of the housing. Remove the cover with a screwdriver. Insert the two 1.5 volt AA batteries in their proper positions. Pay particular attention to polarity. The negative pole is placed against the spring. Equivalent rechargeable batteries (not supplied) may be substituted for the AA's provided. With the AA's one might anticipate approximately 500 hours of service. To enlarge the life time, disable the display backlight.

## 19 Operating of the Contiburette

	<p><b>Attention:</b>  <b>Extreme caution should be exercised when handling toxic, corrosive, fuming, volatile or any other potentially dangerous substances.</b>  <b>Please comply with all safety and accident-prevention regulations applicable to laboratory work.</b></p>
	<p><b>Note:</b>          This instruction sheet does not purport to address all of the safety problems which might result from the use of this instrument, chemicals, reagents, apparatus or equipment employed in any specific test or protocols. It is the responsibility of the user to consult and establish appropriate safety and health practices, and then determine the applicability of regulatory limitations prior to use.          The user himself must determine as to whether the instrument is suitable for the desired purpose. In case of questions, always consult the distributor or manufacturer</p>
	<p><b>Note:</b>          Read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked.</p>
	<p><b>Danger:</b>          Do not use the device in potentially explosive areas. Mortal danger!</p>
	<p><b>Attention:</b>          Use great caution when working in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof. Do not use flammable or explosive substances near the instrument.</p>
	<p><b>Attention:</b>          Please comply with all safety and accident-prevention regulations applicable to laboratory work.</p>
	<p><b>Attention:</b>          Only instructed users may operate the instrument.</p>
	<p><b>Attention:</b>          Observe all markings on the reagent bottles. Dangerous and fuming chemicals must be dispensed in a fume hood. Only employ the instrument for the purpose intended by the manufacturer, and particularly within the resistance limits of the instrument. If in doubt, contact your supplier, or the manufacturer's factory representative at the phone number shown at the front page of this operating instruction.</p>

	<p><b>Attention:</b> Always use the instrument in such a manner that neither the operator, nor any other person is endangered. When dispensing, maintain a physical distance between the instrument and the body. Avoid splashes; use proper connecting vessels, protective clothing and gloves.</p>
	<p><b>Attention:</b> Always check the instrument for leaks and air bubbles. Special attention should be directed to determine that all push-ons, threaded connections and suction tubes are firmly in place before beginning operation. Leaking solutions may endanger persons and materials</p>
	<p><b>Note:</b> In case of trouble (e.g., piston difficult to move or leakage), immediately stop pumping. Clean the instrument according to the cleaning instructions before any further use of the instrument or contact the manufacturer.</p>
	<p><b>Note:</b> Make sure that the unit is only operated under the specified ambient conditions listed under „Technical Data“. Please study the instruction manual before operating the instrument. Please observe all instructions listed in the instruction manual.</p>
	<p><b>Note:</b> Please store the instruction manual in a place easily accessible to every user.</p>
	<p><b>Note:</b> Only use original spares and accessories. These are in alignment with the instrument and ensure reliability.</p>
	<p><b>Note:</b> Please make sure that all safety signs, also on the instrument, are visible.</p>
	<p><b>Attention:</b> Do not open the instrument. Repairs are only to be carried out by trained service technicians.</p>

## 19.1 Switch on/off

The push button to turn on the instrument is located at the front of the burette in the far right position of the group of three. After pushing the button “ON”, the LCD display will be activated. To switch off the instrument, press the button “OFF”.

## 19.2 Filling the instrument and removing the air

After the filling/suction tube has been placed in the liquid, place a beaker or flask under the discharge tube and rotate the hand wheel clockwise until the liquid in the tube is clear of air bubbles. Never turn the hand wheel counter-clockwise.

## 19.3 Operating modes of the Contiburette

Die Contibürette has 2 different operating modes and a set-up menu.

To switch over from one menu to another use the buttons **Left / -** and **Right / +** in the main menu. The actual menu is shown in the menu bar with a black background. To enter a mode/menu or to change a value use the button **Set Value / OK**. This button serves also to confirm a modified parameter. To leave a menu and to refuse a changed value use the **Back / Reset** button.

### 19.3.1 Titrate-Mode

Procedure to titrate:

- Switch on the contiburette.
- Use the buttons **Left / -** and **Right / +** to select the mode **Titrate**.
- The dosed volume is displayed after **Volume:**
- Turning the handwheel effects dosing of the liquid out of the bottle
- The button **Back / Reset** zeroise the dosed volume
- The button **Set Value / OK** recalls a accidentally cleared value.
- The contiburette is switched off with the button **Off**.

### 19.3.2 Dosing-Modus

Additionally to the functions of the Titrate mode, a target volume and the remaining volume to target volume is displayed.

- Switch on the contiburette.
- Use the buttons **Left / -** and **Right / +** to select the mode **Dose**.
- With **Set Value / OK** a target volume is set. The colon after **Dose** changes to an arrow and with **Left / -** and **Right / +** the value is changed. Another **Set Value / OK** confirms the value.
- **Rest:** shows the remaining volume up to the target volume.
- When the target volume is dosed (**Rest: 0.00ml**) , the button **Back / Reset** starts a new cycle, the value **Rest:** changes to the target volume. Another **Back / Reset** zeroes the dosed volume

#### **Attention:**

Even when the target volume is reached, liquid is dispensed when turning the handwheel. The values of target volume and rest volume are only a guidance for the user.

### 19.4 Setup Menu

This menu serves to set:

- Calibration / Adjusting (see Chapt. 20)
- Display backlight (activating the backlight means shorter life time of the batteries)
- Automatic switching off

#### 19.4.1 Setting of the display backlight

The Contiburette has the option to activate a display backlight for better clearness of the display. After a certain time of nonuse the backlight is switched off automatically to spare battery capacity.

Procedure to change the switch off time of the backlight:

- Switch **On** the Contiburette.
- Select **Setup** with the **Left / -** and **Right / +** buttons.
- Enter the submenu with **Set Value / OK** .
- Select with the **Left / -** and **Right / +** buttons **Select> Backlight Time:** and confirm with **Set Value / OK**

- Select with the buttons **Left / -** and **Right / +** the desired time and confirm with **Set Value / OK** .
- To disable the backlight select the time 0 sec. The display shows then **Off**
- Return to the main menu with **Back / Reset**.

#### 19.4.2 Setting of the automatic switching off

The Contiburette has the option to activate a automatically switch off . After a certain time of nonuse the Contiburette is switched off automatically to spare battery capacity.

Proceduere to change the switch off time of the Contiburette:

- Switch **On** the Contiburette.
- Select **Setup** with the **Left / -** and **Right / +** buttons.
- Enter the submenu with **Set Value / OK** .
- Select with the **Left / -** and **Right / +** buttons **Select> Shutoff Time:** and confirm with **Set Value / OK**
- Select with the buttons **Left / -** and **Right / +** the desired time and confirm with **Set Value / OK** .
- To disable the automatically shut off select the time 0 sec. The display shows then **Off**
- Return to the main menu with **Back / Reset**.

## 20 Calibration of the Contiburette

In line with ISO 9000 "Monitoring of Testing Apparatus" and GLP or when using a medium with different density and viscosity other than distilled water there is a possibility to calibrate the Contiburette  $\mu 10$  H. For adjustment of the unit please see chapter 21.

### 20.1 Procedure of calibration:

1. Fill the Contiburette with demi-water, dispense the water into a separate vessel until there are any bubbles are in the dispensing tube left.
2. Set the display to zero.
3. Dispense 5 ml into a vessel.
4. Weigh the dispensed quantity with a precision balance.
5. Calculate the volume, taking the temperature into account.
6. Repeat step 2-5 at least 10 times.
7. Calculate the accuracy A% and coefficient of variation CV% by means of the formulas of the statistical computation.

#### Calculations:

Mean value  $\bar{m} = \frac{\sum m_i}{n}$   $m_i$  : results of weighing,  $n$  : number of weighing

Mean Volume  $\bar{V} = \bar{m} \cdot Z$   $Z$  : Correction factor

Accuracy  $A \% = \frac{\bar{V} - V_0}{V_0}$   $V_0$  : Nominal volume

Coefficient of variation  $CV \% = \frac{100 s}{\bar{V}}$   $s$  : Standard deviation of the results of weighing  $m_i$

A detailed description of this test procedure you find e.g. in DIN EN ISO 8655-6.

Table 1: Correction factors Z ( $\mu\text{l}/\text{mg}$  at 1013 hPa, abstract of EN ISO 8655-6)

Temperature °C	Correction factor Z	Temperature °C	Correction factor Z
15,0	1,0020	23,0	1,0035
15,5	1,0020	23,5	1,0036
16,0	1,0021	24,0	1,0038
16,5	1,0022	24,5	1,0039
17,0	1,0023	25,0	1,0040
17,5	1,0024	25,5	1,0041
18,0	1,0025	26,0	1,0043
18,5	1,0026	26,5	1,0044
19,0	1,0027	27,0	1,0045
19,5	1,0028	27,5	1,0047
20,0	1,0029	28,0	1,0048
20,5	1,0030	28,5	1,0050
21,0	1,0031	29,0	1,0051
21,5	1,0032	29,5	1,0052
22,0	1,0033	30,0	1,0054
22,5	1,0034		

## 21 Adjusting of the Contiburette

The Contiburette is factory adjusted before delivery with water. It might be appropriate to adjust the Contiburette after a longer period of using or for specific applications, e.g. when used ar liquid, which has a significant other density than water.

### 21.1 Procedure to adjust the Contiburette (usercal.):

- Switch **On** the Contiburette.
- Select **Setup** with the **Left / -** and **Right / +** buttons and confirm with **Set Value / OK**.
- Select with the **Left / -** and **Right / +** buttons **Select> Calibration:Factory** and confirm with **Set Value / OK**
- Select **Change Usercal. Yes** with the button **Right / +** and confirm with **Set Value / OK**
- Fill the Contiburette with degased distilled water by turning the handwheel till the tubing is completely free of gas bubbles.
- Select **Fill tubes Yes** with **Left / -** and **Right / +** and confirm with **Set Value / OK**.
- Tare the balance.
- Select **Tare balance Yes** with **Left / -** and **Right / +** and confirm with **Set Value / OK**.

- Turn the handwheel and dose at least the displayed volume. Confirm with **Set Value / OK**.
- Weigh the dosed amount and enter the balance value in grams, using the **Left / -** and **Right / +** buttons. Confirm with **Set Value / OK**.
- Enter the water temperature by using the + and – buttons. Confirm with **Set Value / OK**. The instrument is now doing the automatic calibration itself. The display shows the smallest volume amount you can dose.
- The sign **cal!** is an indication for the user calibration.
- Return to the main menu with **Set Value / OK**.
- With **Back / Reset** you can step backward at any time

## 21.2 Return to Factory Setting

To return to factory setting please carry out the following steps:

- Switch **On** the Contiburette.
- Select **Setup** with the **Left / -** and **Right / +** buttons and confirm with **Set Value / OK**.
- Select with the **Left / -** and **Right / +** buttons **Select> Calibration:User** and confirm with **Set Value / OK**
- Select with the **Left / -** and **Right / +** buttons **Change usercal > Reset** and confirm with **Set Value / OK**
- The sign **cal!** as an indication for the user calibration disappears and the factory calibration is restored.

## 22 Cleaning and Maintenance

This instrument must be cleaned as follows to assure proper functioning and continued accuracy:

- Immediately if the hand wheel becomes sticky or jammed
- Daily after use when these liquids are dispensed:

Solutions prone to crystallization

Alkaline solutions

Organic solvents such as aromatics, chlorinated Hydrocarbons and scintillation liquids

Inorganic oxidising solutions such as buret reagent

- Periodically in order to prolong the life of the instrument.
- When changing the reagent
- Prior to long-term storage

**Caution:** The ceramic parts are subject to binding or freezing if stored after improper cleaning.

**Attention:**

Be careful to avoid any personal injury from used chemicals. While and even after dispensing liquids, the instrument, the filling and the discharge tubes contain the used reagent. Make sure, that during cleaning and maintenance you avoid splashing chemicals. Wear face screens, protective gloves and protective clothes.

### 22.1 Preparation for cleaning

Place the mounted instrument with the bottle into a sink designated for that purpose.

Unscrew the instrument from the bottle and lift the instrument up far enough so that the filling/suction tube is no longer immersed.

Carefully tap the filling tube against the inside of the bottle so that the reagent moves back in. Carefully lift the instrument out of the bottle.

Hold the discharge tube over any designated dispensing receptacle and dispense the remaining reagent.

## 22.2 Cleaning procedure

	<b>Attention:</b> Never turn the hand wheel backward or with force. Be careful to avoid any reagent spillage.
	<b>Attention:</b> Please comply with all safety and accident-prevention regulations applicable to laboratory work.
	<b>Note:</b> The Contiburette $\mu$ 10 C is a measuring instrument and designed to provide high accuracy. To maintain this accuracy, we recommend that this instrument be tested at regular intervals, especially after any mishandling (such as hitting or dropping) of the instrument. Testing of the instrument is provided by the manufacturer for a small fee. Under § 4 of the Weights and Measuring Standards of 12.08.88 (Germany), it is required that regular tests and inspections be performed if the Contiburette $\mu$ 10 C is used as a medical measuring instrument.

## 22.3 Cleaning the Contiburette

- 1 Hold the discharge tube over any designated dispensing receptacle and dispense the remaining reagent (press at the same time "*Manual*" and "+").
- 2 Put the suction tube into cleaning solution designated for that purpose
- 3 Clean the instrument by pumping. We recommend to pump a minimum volume of 50 times the amount of the strokevolume through the pump for a good cleaning.
- 4 Insert the suction tube into distilled water (or other liquids for sterilisation) for rinsing.

type	strokevolume	min. cleaning volume
Contiburette $\mu$ 10 H	200 $\mu$ l	10 ml = 50 turns
Contiburette $\mu$ 20 H	400 $\mu$ l	20 ml = 50 turns

## 22.4 Sterilization

Sterilization is only possible on a chemical basis, such as with the use of alcohol.

## 23 Dismantling and Disposal

### 23.1 Dismantling

**Attention:**

Pumphead and tubing may contain reagents, which endanger persons and material. Make sure of cleaning pumphead and tubing according to Chapter 7 before removing tubing.

1. Switch the instrument off.
2. Disconnect the instrument from the mains.
3. Disconnect the tubing
4. Now the instrument may be removed from the working area.

### 23.2 Disposal

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

## 24 Transport and Storage

### 24.1 Transport/Storage

**Attention:**

Pumphead and tubing may contain reagents, which endanger persons and material. Make sure of cleaning pumphead and tubing according to Chapter 7 before removing tubing.

Prior to transport:

Switch the instrument off and unplug the power supply.

Remove tubing and cables

Do not subject the instrument to mechanical shocks or vibration during transporting.

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)

Store the instrument in a dry environment. Please observe the specified conditions of the ambient:

Ambient temperature: 5-40°C  
Max. relative air humidity: 80%

## 24.2 Return for repair or calibration

	<p><b>Note:</b> For a maximum of protection from health hazards caused by contaminated instruments clean and decontaminate the instrument carefully before returning.</p>
---	---

We intend to give our staff a maximum of protection from health hazards caused by contaminated instruments. We therefore ask for your understanding that we cannot carry out any calibration / repair unless the

### **Declaration on the Absence of Health Hazards**

is submitted completed and signed.

Please copy the declaration in the appendix and attach it completed and signed to the instrument when returned to your distributor or to the manufacturer.

Please provide us with the following supplementary information:

- Detected defect
  
- Media which the instrument has been used with

## 25 Warranty and Liability

The manufacturer agrees to correct for the original user of this product, either by repair, or at the manufacturer's discretion, by replacement, any defects in material or workmanship which develop within 24 months after delivery of this product to the original user. In the event of replacement, the replacement unit will be warranted for the remainder of the original twelve (24) months period of ninety (90) days, whichever is longer.

If this product should require service, contact your distributor or manufacturer for necessary instructions.

This warranty shall not apply if the defect or malfunction was caused by accident, neglect, unreasonable use, improper service, or other causes not arising out of defects in material or workmanship.

There are no warranties, expressed or implied, including, but not limited to, those of merchantability or fitness for a particular purpose, which extended beyond the description and period set forth herein.

The manufacturer's sole obligation under this warranty is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind resulting from improper use or misuse or possession of the product.

	<p><b>Attention:</b> In case of malfunction do not try to carry out any repair works. The instrument does not consist of any part which may be serviced or maintained by the user. Any attempt by the user to repair the unit will cancel the warranty.</p>
	<p><b>Attention:</b> Do not open the instrument. Any work on the electronics of the unit should only be done by knowledgeable and trained personnel.</p>

## 26 Technical Data

**Attention:**

The user has to determine, if the instrument is suitable for his specific application. If there are any further questions, contact your local dealer or the manufacturer.

EX 20 °C

Subdivision	0,01 ml ( $\mu$ 10 H) 0,02 ml ( $\mu$ 20 H)
Accuracy A	$\leq$ 0,2%
Coefficient of variation CV	$\leq$ 0,1%
Error limits according to DIN ISO 8665-3 relative to the nominal volume or $\pm$ 10 $\mu$ l	
Operating Temperature Range:	10 – 40°C
Transport and Storing Temps:	Maximum 60°C
Power supply:	2 Batteries type AA (LR6, AM3)

### 26.1 Chemical Resistance:

Materials which come into contact with the dispensed Liquids:

**Aluminium oxide Al<sub>2</sub>O<sub>3</sub> (99,7%)**

Sapphire hard, fused crystalline alumina - excellent chemical resistance to most commonly used fluids - mechanically resistant to common abrasives.

**ETFE & FEP**

Excellent chemical resistance characteristics.

## 27 Declaration on the Absence of Health Hazards

Please copy this declaration and attach it completed and signed to the instrument

Device designation:.....

Serial No.: .....

### The Undersigned hereby declares:

- ◆ That the instruments have been carefully cleaned and decontaminated before shipment.
- ◆ That the instruments pose no danger through bacteriological, chemical, radiological or viral contamination.
- ◆ To be authorised to make declarations on behalf of the Institution represented.
- ◆ That he / she is aware that shipment of contaminated instruments is a violation of law, and that he / she personally and the Institution represented may be held liable for any damages caused by contaminated instruments.
- ◆ For calibrating service only: minor repairs of a value up to € 30,--+ VAT will be carried out and invoiced without further queries (cross out if not applicable).

Sender:

Firm / Laboratory:

Name

.....

.....

Position

Address:

.....

Date, Signature

.....

Tel. for enquiry:

.....

.....

- ◆ In case of Return for Repair, please provide us with the following supplementary information:

Detected defect:

.....

Media which the instrument has been used with:

.....