**CAT**

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
*T6, T10, T17, T/G20, T/G30, T/G40*

## Content

1	General Information .....	4
1.1	Combinations of Generators (using the example of shaft T20) .....	4
2	Unpacking the Instrument .....	5
2.1	Scope of delivery .....	5
2.2	Homogenising Tools and Accessories (to be ordered separately): .....	6
3	Safety Information .....	7
3.1	Explanation of the safety information in the instruction manual .....	7
3.2	General Safety Information .....	7
4	Description of a Homogenising Tool .....	8
4.1	Exchanging the Generator: .....	8
4.1.1	Disassembling Shaft T6.....	9
4.1.2	Reassembling Shaft T6.....	10
4.1.3	Disassembling Shaft T10/T17 .....	10
4.1.4	Reassembling Shaft T10/T17 .....	11
4.1.5	Disassembling Shaft T20/T30/T40 .....	12
4.1.6	Reassembling Shaft T20/T30/T40 .....	13
4.1.7	Disassembling Shaft G20/G30 .....	13
4.1.8	Reassembly of Shaft G20/G30 .....	14
4.1.9	Disassembling Shaft G40.....	15
5	Maintenance and Cleaning .....	16
5.1	Cleaning the homogenising tool.....	16
5.2	Maintenance homogenizing tools.....	16
6	Transport and Storage.....	16
7	Disposal .....	17
8	Warranty and Liability.....	17
9	Technical Data.....	18
9.1	Technical Data.....	18
9.2	Technical Data Rotors and Stators.....	19
10	Repairs .....	19
10.1	Repair Return Form.....	20
1	Allgemeine Information .....	22
1.1	Generator Kombinationen .....	22
2	Auspacken des Gerätes .....	23
2.1	Lieferumfang und Zubehör.....	23
2.2	Dispergierwerkzeuge und Zubehör (separat zu bestellen):.....	24

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
## 1 General Information

A dispersion tool consists of a shaft and a generator (rotor and stator).

### Shaft:

The shaft mainly consists of the shaft tube, the drive shaft rotating shaft inside the shaft tube, the bearing and the shaft sealing.


The shaft cannot be used alone. Only combinations of shaft and generator (A) can be used for working with a medium.


### Generator:

The generator (A) consists of a rotor (1) and a stator (2). These parts can be unscrewed from the shaft.


#### 1.1 Combinations of Generators (using the example of shaft T20)

To create an F-type generator use an F-rotor plus an F-stator. This combination is suitable for aqueous media.

To create an N-type generator, use an F-rotor plus a V-stator. This combination is suitable for normal media.

To create a V-type generator use a V-rotor plus a V-stator. This combination is suitable for coarse media.

**Note:**

*Rotors and stators form matching pairs. They can be combined to form F-, N-, V-versions. Exceptions are the knife generators. M-rotors and M-stators form pairs and must not be confused.*

## 2 Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport to be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

### 2.1 Scope of delivery

The dispersion tools you have ordered are always supplied fully assembled.

Please check that the package contains the following:

- 1 fully assembled Shaft
- 1 Instruction Manual

## 2.2 Homogenising Tools and Accessories (to be ordered separately):

Tools	Description	Part No.
T 6	Shaft, 6 mm diameter, 55 mm long, complete with generator	60420*
T 6	Shaft, 6 mm diameter, 120 mm long, complete with generator	60410*
T 10	Shaft, 10 mm diameter complete with generator	60421*
T 17	Shaft, 10 mm diameter, complete with 17 mm generator	60426*
T 20	Shaft, 20 mm diameter complete with generator	60422*
G 20	Shaft, 20 mm diameter complete with generator	60423*
T 30	Shaft, 20 mm diameter complete with 30 mm generator	60424*
G 30	Shaft, 20 mm diameter complete with 30 mm generator	60425*
T 40	Shaft, 20 mm diameter complete with 40 mm generator	60414*
G 40	Shaft, 20 mm diameter complete with 40 mm generator	60415*
Universal key	Universal rotor wrench for Shafts 6 - 30 mm	60470-00
Socket wrench	for all shafts	60471-00
Rotor wrench	for shaft T 10	60472-00
Rotor wrench	for shaft T 20/G20	60473-00
Rotor wrench	for shaft T 30/G30	60476-00
Rotor wrench	for shaft T 40	60477-00

### **Important Note:**

*Read this instruction manual carefully before operating the instrument. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer at the following address*

# **CAT**

**Ingenieurbüro CAT  
M. Zipperer GmbH  
Etzenbach 16  
D-79219 Staufen  
+49-7636-7803-0**

### **Important Note:**

*Please store the instruction manual in a place easily accessible to every user.*

### **Important Note:**

*Make sure that the unit is to be operated under the conditions corresponding to that as described under the section titled "Technical Data"*

### 3 Safety Information

#### 3.1 Explanation of the safety information in the instruction manual

*This symbol emphasizes the following notices, which have to be strictly observed by the user! Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.*

*This symbol emphasizes notices, which have to be observed by the user to assure a safe operation of the unit.*

#### 3.2 General Safety Information

**Attention:**

*Please comply with all safety and accident-prevention regulations applicable to laboratory work.*

**Attention:**

*Only instructed users may operate the instrument.*

**Important Note:**

*Make sure that all safety information of the instrument respectively on the device itself are clearly visible during operation*

**Attention:**

*If noise level exceeds 85 dB (A) at the work station operator must use ear plug and/or ear mufflers.*


*The Generators (A) have sharp edges. Handle with care. Risk of injuries!*

## 4 Description of a Homogenising Tool

The homogenizing tools come always assembled. The connecting adapter (drive pin (4)) is located on the top of the shaft and is used to attach the shaft to the drive motor or to the socket wrench (B) during assembling or disassembling.

On the other end of the shaft there are stator and rotor. Socket wrench and rotor wrench (C) are used for assembly resp. disassembly.


**Never run homogenising tools dry** as the bearings and gaskets will be damaged if the generator is not cooled by the medium. The distance between the homogenising tool and the vessel bottom should not be less than 10mm. The filling of the medium may not be less than about 55 mm. To avoid a vortex it is necessary to insert the dispersing shaft out of the centre of the vessel. If these conditions are met the unit is ready for operation.

### 4.1 Exchanging the Generator:


**Attention:**

*The generators (A) have sharp edges. Handle with care. Risk of injuries!*

Socket Wrench (B)

Universal Rotor Wrench (C)

For undoing the rotor (1), use the two special keys (to be ordered separately). The socket wrench (B) is used to counter-hold the Shaft, and the rotor wrench (C) to undo the rotor (1).

To exchange the stator (2), remember that the fastening thread is **counter-clockwise**.


**Important Note**

*Always unscrew the rotor (1) first before unscrewing the stator (2).*


### 4.1.1 Disassembling Shaft T6



#### Shaft T6

Bearing at the bottom of the shaft  
Lager am Schaffende

Bearing at the top of the shaft  
Lager antriebsseitig

- 01 Rotor
- 02 Stator
- 03 Rotating Shaft
- 04 Drive Pin
- 05 Threaded Ring
- 06 PTFE Bearing

- 07 Shaft Tube Bearing
- 08 Thread Pin
- 09 Ball Bearing
- 10 PTFE Disc

The shaft comes always fully assembled. The socket wrench (B) is used to counter hold the rotating shaft (3) and to undo the nut (D). The rotor wrench (C) is used to loosen respectively to fasten the rotor (1). To loosen the stator shaft you need a 2.5 Allen key.

To disassemble the generator (A), carefully insert the rotor wrench (C) into the rotor (1). Now place the drive pin (4) located at the rotating shaft (3) into the socket wrench (B). Carefully turn to the left to unscrew the rotor (1).

The removal of the bearing is only necessary for cleaning or replacement purposes.

#### 4.1.2 Reassembling Shaft T6

Before assembling, be sure the bearing is in place. Insert the socket wrench (B) on top of the shaft as explained above. Hold the shaft and socket wrench in one hand and screw the rotor (1) on clockwise (to the right). Tighten the rotor (1) using the rotor wrench (B).

Insert the shaft into the drive motor firmly with a "click". Tighten the shaft securely with the screw on the motor flansk.


#### **Important**

*NEVER change the rotor /stator while the shaft is attached to the motor! NEVER run the motor dry. ALWAYS insert the shaft into the media BEFORE turning on the motor!*

#### 4.1.3 Disassembling Shaft T10/T17

##### Shaft T10/T17

6 PTFE Bearing  
7 Sliding Sleeve  
8 Spring  
3 Rotating Shaft  
4 Drive Pin

11 Shaft Tube  
9 Ball Bearing  
5 Threaded Ring  
10 PTFE Disc

The shaft comes always fully assembled. The socket wrench (B) is used to counter hold the rotating shaft (3) and to undo the nut (I). The rotor wrench (C) is used to loosen respectively to fasten the rotor (1).

To disassemble the generator (A), carefully insert the rotor wrench (C) into the rotor (1). Now place the drive pin (4) located at the rotating shaft (3) into the socket wrench (B). Carefully turn to the left to unscrew the rotor (1). Take care not to bend the teeth of the rotor (1). The stator (2) can now be unscrewed turning to the right by hand.

***Please note***

*Left hand thread.*

The PTFE bearing (6) can now be pulled out. The bearing section consists of two parts, the bearing itself and the shaft sleeve.

The removal of the bearing is only necessary for cleaning or replacement purposes.

#### **4.1.4 Reassembling Shaft T10/T17**

Before assembling make sure the bearing is in place. Screw on the stator (2) in counter clockwise direction. Insert the socket wrench (B) on top of the shaft as described under "disassembling". Hold the shaft and socket wrench (B) in one hand and screw on the rotor (1) (clockwise -to the right). Tighten the rotor (1) using the rotor wrench (C).

Insert the shaft into the drive motor firmly with a "click". Tighten the shaft securely with the screw on the motor flansk.


***Important***

*NEVER change the rotor /stator while the shaft is attached to the motor! NEVER run the motor dry. ALWAYS insert the shaft into the media BEFORE turning on the motor!*

## 4.1.5 Disassembling Shaft T20/T30/T40

### Shaft T20/T30/T40

- |                    |                   |
|--------------------|-------------------|
| 1 Rotor            | 3 Rotating Shaft  |
| 2 Stator           | 12 Retaining Disc |
| 7 Sliding Sleeve   | 10 PTFE Disc      |
| 14 Bearing Bushing | 9 Ball Bearing    |
| 6 PTFE Bearing     | 5 Threaded Ring   |
| 13 Retaining Disc  |                   |
| 11 Shaft Tube      |                   |

#### 4.1.5.1 Disassembly of bearing at the bottom of the shaft

The homogenizer tool comes always assembled. To disassemble hold the shaft in one hand and unscrew the stator (2) clockwise (to the right). Insert the socket wrench (B) at the top of the shaft by placing the drive pin (4) into the large slot provided by the socket wrench. Turn the socket wrench (B) until the two small rods at the end of the socket wrench (B) go into the two small holes located at the top of the shaft. With one hand holding the shaft and the socket wrench take the rotor wrench (C) in the other hand and place the flat bar between the teeth of the rotor (1). Turn counter-clockwise (to the left). The rotor (1) now can be unscrewed easily. Now unscrew the stator (2) clockwise (to the right). The bearing now can be pulled out.

The removal of the bearing is only necessary for cleaning or replacement purposes.

The bearing assembly at the bottom of the shaft consists of four parts: Sliding sleeve (7), bearing bushing (14), PTFE bearing (6) and a retaining disc (13).

#### 4.1.5.2 Disassembly of the bearing at the top of the shaft

The bearing assembly at the top of the shaft consists of four parts. Retaining disc (12), PTFE disc (10), ball bearing (9) and threaded ring (5).

Reassemble the bearing in reverse sequence and securely fasten the threaded ring (5).

#### 4.1.6 Reassembling Shaft T20/T30/T40

Before assembling make sure the bearing is in place. Screw on the stator (2) in counter clockwise direction. Insert the socket wrench (B) on top of the shaft as described under "disassembling". Hold the shaft and socket wrench (B) in one hand and screw on the rotor (1) (clockwise = to the right). Tighten the rotor (1) using the rotor wrench (C).

Insert the shaft into the drive motor firmly with a "click". Tighten the shaft securely with the screw on the motor flange.


#### **Important**

*NEVER change the rotor /stator while the shaft is attached to the motor! NEVER run the motor dry. ALWAYS insert the shaft into the media BEFORE turning on the motor!*

#### 4.1.7 Disassembling Shaft G20/G30

##### Shaft G20/G30

- 16 Spiral Spring
- 15 Retaining Disc
- 14 O-Ring
- 13 Slip Ring
- 12 Slip Disc
- 7 O-Ring (x2)
- 10 Spacer Disc
- 6 PTFE Disc (x2)
- 9 Ball Bearing

- 3 Rotating Shaft 3-part
- 4 Drive Pin (not shown in picture)
- 11 Shaft Tube
- 9 Ball Bearing
- 8 Retaining Ring
- 5 Threaded Ring

#### 4.1.7.1 Disassembly of bearing at the bottom of the shaft G20/G30

The homogenizer tool comes always assembled. To disassemble hold the shaft in one hand and unscrew the stator (2) clockwise (to the right). Insert the socket wrench (B) at the top of the shaft by placing the drive pin (4) into the large slot provided by the socket wrench (B). Turn the socket wrench (B) until the two small rods at the end of the socket wrench (B) go into the two small holes located at the top of the shaft. With one hand holding the shaft and the socket wrench (B) take the rotor wrench (C) in the other hand and place the flat bar between the teeth of the rotor (1). Turn counter-clockwise (to the left). The rotor (1) now can be unscrewed easily. Now unscrew the stator (2) clockwise (to the right). The bearing now can be pulled out. The removal of the bearing is only necessary for cleaning or replacement purposes.

The bearing assembly at the bottom of the shaft consists of eleven Parts: spiral spring (16), retaining disc (15), O-ring (14), ceramic slip ring (13), slip disc (12), O-ring (7), spacer disc (10), two PTFE discs (6) and a ball bearing (9).

#### 4.1.7.2 Disassembly of the bearing at the top of the shaft

The bearing assembly at the top of the shaft consists of three parts. Ball bearing (9), retaining ring (8) and threaded ring (5).

Reassemble the bearing in reverse sequence and securely fasten the threaded ring (5).

The unit comes always assembled. The socket wrench (B) is used to counter-hold the drive shaft (3). The rotor wrench (C) is used to unscrew respectively to fasten the rotor (1).

To disassemble the generator (A), carefully insert the rotor wrench (C) into the rotor (1) and place the drive pin (4) located at the drive shaft (3) into the socket wrench (B). Now carefully turn to the left to unscrew the rotor (1). Take care not to bend the teeth of the rotor (1).

The removal of the bearing should only be necessary for cleaning or replacement purposes.

#### 4.1.8 Reassembly of Shaft G20/G30

Before assembling make sure the bearing is in place. Screw on the stator (2) in counter clockwise direction. Insert the socket wrench (B) on top of the shaft as described under "disassembling". Hold the shaft and socket wrench (B) in one hand and screw on the rotor (1) (clockwise -to the right). Tighten the rotor (1) using the rotor wrench (C).

Insert the shaft into the drive motor firmly with a "click". Tighten the shaft securely with the screw on the motor flansk.


#### **Important**

*NEVER change the rotor /stator while the shaft is attached to the motor! NEVER run the motor dry. ALWAYS insert the shaft into the media BEFORE turning on the motor!*

## 4.1.9 Disassembling Shaft G40

### Shaft G40

03 Rotating Shaft (three part)

04 Drive Pin

05 Threaded Ring

06 Plastic Gasket

07 O-Ring

08 Retaining Ring

09 Ball Bearing

10 Spacer Disc

11 Shaft Tube

12 Ceramic Slip Disc

13 Ceramic Slip Ring

14 Slip Ring Disc

15 O-Ring

16 Spiral Spring

### 4.1.9.1 Disassembly of bearing at the bottom of the shaft G20/G30

The homogenizer tool comes always assembled. To disassemble hold the shaft in one hand and unscrew the stator (2) clockwise (to the right). Insert the socket wrench (B) at the top of the shaft by placing the drive pin (4) into the large slot provided by the socket wrench (B). Turn the socket wrench (B) until the two small rods at the end of the socket wrench (B) go into the two small holes located at the top of the shaft. With one hand holding the shaft and the socket wrench (B) take the rotor wrench (C) in the other hand and place the flat bar between the teeth of the rotor (1). Turn counter-clockwise (to the left). The rotor (1) now can be unscrewed easily. Now unscrew the stator (2) clockwise (to the right). The bearing now can be pulled out. The removal of the bearing is only necessary for cleaning or replacement purposes.

The bearing assembly at the bottom of the shaft consists of eleven Parts: spiral spring (16), retaining disc (15), O-ring (14), ceramic slip ring (13), ceramic slip disc (12), two O-rings (7), spacer disc (10), two PTFE discs (6) and a ball bearing (9).

## 5 Maintenance and Cleaning

### 5.1 Cleaning the homogenising tool

**Attention:**

*The Generators have sharp edges. Handle with care. Risk of injuries!*

To avoid clogging of the shaft, clean shaft, generator and in case of G-shaft the sealing parts after each use. This is done by operating it in a solvent which dissolves substance residues and is not harmful to the gasket. This is usually sufficient to clean the Generator.

**Chemical sterilization** may be also a method. General-purpose disinfectants such as formalin, alcohol, etc. may be used. It is important to remove disinfectant residues with sterilized water.

**Important Note:**

Make sure that the bearings, O-rings and gaskets are resistant to solvents.

**Sterilization with moist heat:**

This denotes the use of a steam jet pressurized to 2 bars at 120°C.

### 5.2 Maintenance homogenizing tools

The gaskets in the homogenizing tools must be constantly monitored. In the event of leakage the suction effect of the rotating Shaft can cause the medium to penetrate as far as the drive unit. If liquid emerges from the side hole at the top of the Shaft tube stop work immediately and check the gaskets. The function of the Generators depends on the condition of the sharp edges on the rotor and stator. These edges may be blunted very quickly in abrasive media, reducing the effectiveness of homogenizing.

## 6 Transport and Storage

**Prior to transport:**

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)
- Storage ambient:

Max. ambient temperature : RT to +40°C

Max. humidity: 80%

## 7 Disposal

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: Metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

## 8 Warranty and Liability

The manufacturer agrees to either repair, or replace, at the manufacturer's discretion, any defects in materials or workmanship which develop within 24 months of the delivery of this product to the original user. In the event of replacement, the replacement unit will be guaranteed for the remainder of the original twenty-four (24) months period or ninety (90) days, whichever is longer.

If this product should require service, contact your local distributor or manufacturer for necessary instructions.

This guarantee will not apply if the defect or malfunction was caused by accident, neglect, unreasonable use or fitness for a particular purpose, which extend beyond the description and period set forth herein.

The manufacturer's sole obligation under this guarantee is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind, resulting from use or possession of the product.

### **Attention:**

The user has to determine, if the instrument is suitable for his specific application. If there are any further queries, contact your local dealer or the manufacturer direct.

## 9 Technical Data

### 9.1 Technical Data

Shaft Type	Generator $\varnothing$ in [mm]	Shaft Length in [mm] incl. generator	Immersion Depth in [mm]	Volume in [ml]	Materials
T6	6	105	55	0.1 - 50	Stainless Steel, PTFE
T6 long	6	170	120	0.1 - 50	Stainless Steel, PTFE
T10	10	205	150	1.0 - 250	Stainless Steel, PTFE
T17	17	205	150	5.0 - 1,000	Stainless Steel; PTFE
T20	20	265	200	10.0 - 2,000	Stainless Steel, PTFE
G20	20	265	215	10.0 - 2,000	Stainless Steel, PTFE
T30	30	255	190	30.0 - 5,000	Stainless Steel, PTFE
G30	30	255	200	30.0 - 5,000	Stainless Steel, PTFE
T40	40	360	250	100.0 - 20,000	Stainless Steel, PTFE
G40	40	360	280	100.0 - 20,000	Stainless Steel, PTFE

## 9.2 Technical Data Rotors and Stators

Description	Teeth	Outer $\varnothing$ in [mm]	Inner $\varnothing$ in [mm]	Part no.
Stator 6 V	7	6	4.2	60450-00
Stator 6 F	11	6	4.2	60451-00
Stator 6 V	7	6	4.2	60452-00
Stator 6 F	11	6	4.2	60453-00
Rotor 6 V	2	4	-	60460-00
Rotor 6 F	4	4	-	60461-00
Stator 10 V	7	10	8.15	60430-00
Stator 10 F	17	10	8.15	60431-00
Rotor 10 V	2	7.5	-	60440-00
Rotor 10 F	4	7.5	-	60441-00
Stator 17 V	9	17	13.3	60438-00
Stator 17 F	13	17	13.3	60439-00
Rotor 17 V	2	12.5	-	60448-00
Rotor 17 F	6	12.5	-	60449-00
Stator 20 V	9	20	16	60432-00
Stator 20 F	13	20	16	60434-00
Rotor 20 V	2	15	-	60442-00
Rotor 20 F	6	15	-	60444-00
Stator 30 V	13	30	26	60435-00
Stator 30 F	17	30	26	60437-00
Rotor 30 V	3	25	-	60445-00
Rotor 30 F	6	25	-	60447-00
Stator 40 V	13	40	35	60456-00
Stator 40 F	23	40	35	60457-00
Rotor 40 V	4	33	-	60466-00
Rotor 40 F	8	33	-	60467-00

## 10 Repairs

*When returning instruments for repair that have come into contact with hazardous substances, please:*

*Fill in attached "Repair Return Form"*

*Provide precise information on the relevant medium*

*Take protective measures to ensure the safety of our receiving and maintenance personnel. Mark the package as appropriate for hazardous materials.*

## 10.1 Repair Return Form

### CONTACT/USER INFORMATION

Contact:	Phone No.
Fax No.	Email:
Billing:	Shipping:
Company	Company
Address	Address

### INSTRUMENT INFORMATION

Model	Serial No.	
Please describe all problems/malfunctions		
Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

\*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to CAT and an MSDS for hazardous material must be included with the instrument.

### RETURN SHIPPING

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	--	------------------------------------	----------------------------------

\*\*Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

### PACKAGING INSTRUCTIONS TO RETURN A CAT INSTRUMENT FOR REPAIR

- ✓ Remove all accessories (e.g. homogniser tools, stirring paddles) from the instrument
- ✓ Clean excess testing material off the instrument/accessory
- ✓ Include MSDS sheets for all hazardous materials used with this instrument
- ✓ Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- ✓ DO NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the instrument. (Applicable for overhead stirrers and homogenisers.)
- ✓ Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.

Example: "Handle with care" or "Fragile- Delicate Instrument"

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**Ingenieurbüro CAT, M. Zipperer GmbH, Service Department, Etzenbach 16, D-79219 Staufen**